
*Loveline*® & flex

## ENGLISH

Enjoy your *Loveline*<sup>®</sup>  
but first please take a moment  
to review this manual. ;-)

### ♥ FLEX What's included


- ♥ Set of 3 Kegel balls of varying weights and sizes
- ♥ Satin privacy pouch

### ♥ Why FLEX?

FLEX is a graduated Kegel exercise aid. We developed this Kegel weight set to accommodate women of varying sizes and PC muscle strength. What is a PC muscle? The “PC” (or Pubococcygeus) muscle stretches from the inside of the vagina to the spinal cord. The PC muscles help to keep your bladder, uterus and other internal organs in place. They may also assist in the control of opening and closing your urethra, the tube through which urine leaves your body. Kegel exercises consist of repeatedly contracting and relaxing the muscles that form part of the pelvic floor. Kegel exercises on a regular basis may help to tone and strengthen your pelvic floor muscles.

**Note:** Always consult your physician if you are pregnant or have any medical concerns prior to using this product. No medical claims are warranted or implied by the use of this product. For adult use only.


## ♥ Using FLEX?

Inserting FLEX balls is similar to inserting a tampon.

- To ease insertion, use a small amount of water-based lubricant on the FLEX ball you plan to use. Silicone based lubricant is NOT recommended for use with this product as silicone based lubricant can deteriorate the silicone exterior on the FLEX product itself.
- Insert the Kegel ball into your vagina so it rests above your pelvic floor muscle.
- During insertion and positioning of FLEX make sure the retrieval cord is outside of the vagina for easy removal.

To remove FLEX, simply wrap the retrieval cord around your finger and gently tug on the cord to remove from the vagina.

## ♥ Kegel Exercise Ideas with FLEX

### **Step 1 – Discover your pelvic floor muscles.**

If this is your first time using Kegel weights, we recommend that you start with our lightest and easiest Kegel ball. This happens to also be the largest one. Its larger size makes it easier for the user to hold inside. Insert the Kegel ball (as described earlier) then tightly squeeze your pelvic floor muscles while sitting. Then relax and try repeating this a few times while you get acquainted with your pelvic muscles. If you have difficulty finding and contracting your pelvic floor muscle, keep trying until you feel the sensation of lifting and squeezing. Each woman's body is different, so at the beginning, just try feeling your PC muscle tensing and relaxing. Having control is a great start!

### **Step 2 – Develop a daily routine or PC workout program.**

Many women integrate exercise of their PC muscles as part of their daily routine. Some habitually do them each day, for example, on their drive to work, while others exercise multiple times a day, whenever they remember.

### **Step 3 – Start FLEXing**

Try keeping your pelvic muscles tense for 5 seconds while sitting. Then relax for 5 seconds. Holding FLEX for 5 seconds should be your goal, so don't be discouraged if you are unable to hold it that long at first. Although you may be able to do only one Kegel exercise when you first start out, over time your muscles will become stronger and more accustomed to doing more repetitions and carrying a heavier weighted ball. Now, try holding FLEX for 10 seconds following with 10 seconds of relaxation. Many experts recommend at least 3 sets of 10 repetitions at least once daily. When it becomes easier, try increasing your repetitions to 15. There's no pressure as to how often or how many you do! When you start to feel more control, graduate to heavier weights or even to the two-ball set.


## ♥ Specifications

♥ **Material:** Silicone

### Beginner

♥ **Weight:**  
35g (1.2oz)


♥ **Diameter:**  
36 mm (1.4")

### Intermediate

♥ **Weight:**  
45g (1.6oz)

♥ **Diameter:**  
29 mm (1.1")

### Advanced

♥ **Weight:**  
85g (3.0oz)

♥ **Diameter  
(of each heart):**  
29 mm (1.1")

## ♥ Care & Cleaning

- Clean after every use with mild soap and water. Dry thoroughly.
- Store in privacy pouch to keep your Lovelife clean and safe until next use.

## ♥ Important Safeguards

- Stop using if you experience pain or discomfort.
- Keep away from children

## ♥ Warranty

OhMiBod warrants this product for a period of one (1) year, from the date of original purchase, against defects due to faulty workmanship or materials. In the event that this product fails during the warranty period, simply email us at [customercare@ohmibod.com](mailto:customercare@ohmibod.com) and provide us with a valid proof of the date of original purchase and we will replace the faulty product at no extra charge.

This warranty does not cover cosmetic deterioration caused by normal wear and tear or damage caused by accident, misuse, neglect or other use not in accordance with this instruction manual. Any attempt to open or take apart this product or its accessories will void the warranty.

## ♥ Contact Us

If you have any questions regarding your Lovelife product, please check out additional information at [www.ohmibod.com](http://www.ohmibod.com). You can also reach us via email at [customercare@ohmibod.com](mailto:customercare@ohmibod.com).

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Copyright© Suki, LLC 2018. All rights reserved.